

Quantum-Si's Next-Generation Protein Sequencer™ Enables Protein Detection and Peptide Characterization from Biological Samples

P-III-0753

John Kudolo, Khanh D. Q. Nguyen, Kenneth Skinner, Mathivanan Chinnaraj, Abde Ali Kagalwalla, Ben Moree, and Brian Reed

Quantum-Si Incorporated, 29 Business Park Drive, Branford, CT 06405, U.S.A.

INTRODUCTION


Studies of the proteome would greatly benefit from methods that enable direct protein detection and peptide characterization with exceptional sensitivity. In this context, we present the unique capabilities of Next-Generation Protein Sequencing™ on Quantum-Si's Platinum® instrument, which utilizes single-molecule protein sequencing to achieve detection of unknown proteins and protein variants. Individual peptides are digested from proteins, immobilized on a semiconductor chip, and probed by dye-labeled N-terminal amino acid (NAA) recognizers, followed by subsequent aminopeptidase cleavage to expose each NAA in the peptide for recognition. Through the recording and analysis of fluorescent intensity, lifetime, and binding kinetics of each NAA binding event, we can successfully identify proteins and detect PTMs.

To demonstrate the power of protein sequencing, we identified proteins from serum without prior knowledge, identified pathologic PTMs, and utilized protein barcodes to screen and characterize proteins. We isolated proteins from human serum via immunoprecipitation or SDS-PAGE and correctly identified them from the sequencing data with high confidence by mapping to an 8,000-protein reference panel. Additionally, we demonstrated the power of Platinum to detect PTMs based on kinetic changes by detecting citrullination and dimethylation of arginine—two PTMs that play key roles in disease states such as cardiovascular disease, autoimmune disease, and cancer. Finally, we showcased the detection of peptide biomarkers along with the utilization of barcoding techniques to streamline protein engineering, antibody engineering, and enzyme engineering.

These results demonstrate the transformative potential of single-molecule protein sequencing using Platinum for comprehensive proteome analysis. Our technology enables the detection and characterization of proteins and peptides in complex mixtures and biofluids, as well as the detection of critical PTMs implicated in human pathophysiology. As proteomics continues to advance, Quantum-Si's platform offers a promising approach to drive new insights into protein research and disease mechanisms.

METHODS

- Proteins are reduced, alkylated, and digested with LysC.
- Peptides are functionalized, conjugated, and immobilized on the surface of a proprietary semiconductor chip.
- Fluorescently labeled N-terminal amino acid (NAA) recognizers and aminopeptidases are added to the semiconductor chip.
- Fluorescent intensity and duration of each NAA binding event generates a unique kinetic signature.
- Kinetic signatures are analyzed to align reads to reference peptides and compute false discovery rate (FDR).




RESULTS


Unknown Proteins from Human Serum Were Correctly Identified with Platinum

Human serum was run on SDS-PAGE, and individual bands were excised, followed by in-gel digestion, library preparation, and sequencing on Platinum. **Three proteins, alipoprotein A-I, albumin, and transferrin were correctly identified** without prior knowledge using the Protein Inference analysis workflow. The results were validated with Western blot.


Below, the 25-kDa band was excised, sequenced, and identified to be APOA1.




Similarly, the 55-kDa band was identified to be ALBU with 99.99% confidence.



Finally, the 75-kDa band was identified to be TRFE with 99.99% confidence.






Below is a detailed analysis of the TRFE dataset, showing the coverage needed to correctly identify without prior knowledge, together with the number of alignments, level of confidence, and kinetic signatures of each peptide.




Kinetic Signatures Enabled the Detection of Citrullination and Dimethylation of Arginine


Asymmetric dimethylation of arginine (ADMA) led to an increase in pulse duration (PD) of the preceding Y residue. Symmetric dimethylation of arginine (SDMA) abolished arginine recognition by the R recognizer.



Citrullination of arginine abolished arginine recognition by the R recognizer and led to an increase in PD of the preceding L residue.



When spiked into a protein sample, a similar result was observed. Citrullination abolished arginine recognition and led to the recognition of the preceding L residue due to an increase in PD.




Protein Barcodes Were Detected at Sub-Picomolar Concentration

Protein barcodes are short peptides that are highly visible and distinguishable with NGPS on Platinum. When expressed with proteins of interest, they can facilitate multiplexing of functional expression of proteins.

These barcodes can be used for:


- Directly correlate multiple protein functions to sequence at once to increase throughput.
- Study protein trafficking of proteins from various organelles.
- Identify and characterize protein-protein interactions.
- Screen and characterize proteins with different properties and functions (e.g. mRNA vaccines)



Mix of Engineered Proteins



In this experiment, we demonstrated the detection of these barcodes in a mixture of a wide dynamic range from 1,800–0.18 pmol. All barcodes were detected with false discovery rates (FDR) of < 2%.



REFERENCES

Brian D. Reed et al, *Science* 2022, 378 (6166) 186-192.

All trademarks are the property of Quantum-Si, Inc. or their respective owners. For specific trademark information, see www.quantum-si.com/privacy-policy.

Research use only. Not for use in diagnostic procedures.