Single-molecule sequencing of immunopeptides and antigenic protein markers on the Platinum platform Kenneth A. Skinner^{1#}; Saketh Kapoor^{2#}; Étienne Caron² **QUANTUM SI** 1Quantum-Si Incorporated, Branford, CT; ²Department of Immunobiology, Yale University, New Haven, CT #Equal contribution

Abstract

- Motivation: A single peptide-MHC complex can activate immune cells; thus, single-molecule sequencing has the potential to provide the requisite sensitivity to elucidate the immunopeptidome and profile tumor-associated antigens.
- Benchtop platform: Platinum reveals the primary structure of single protein molecules, interrogating immunopeptides and proteins that elicit immune responses.

Platinum workflow:

- Derivatized peptides are conjugated via C-termini and loaded into nanoscale reaction chambers of a semiconductor chip.
- During sequencing, N-terminal amino acid (NAA) recognizers, labelled with different fluorophores, reversibly bind to cognate NAAs and generate distinct binding patterns.
- Aminopeptidases sequentially cleave NAAs to expose subsequent amino acids for recognition. Upon completion of a sequencing run, data is automatically transmitted for **cloud-based analysis**.
- Traces are segmented into contiguous groups of pulses, or recognition segments (RSs), each representing interactions between recognizers and NAAs.
- The order of recognizer binding and kinetic properties of RSs are analyzed to generate highconfidence alignments to reference peptides and identify proteins.
- **Results:** Platinum detects peptides selected from a mixture of viral epitopes that activate T cells.
 - Platinum not only sequences cytokines that can remodel the immunopeptidome but also reveal the constituents of immunoprotein mixtures.
 - Platinum enables the detection of sites subjected to **post-translational modification (PTM)**, such as citrullination of arginine residues.
- **Summary:** Platinum profiles the amino acid sequence of viral peptides and proteins integral to immune checkpoint.

Next-Generation Protein Sequencing[™]



PD (s) 0.45 2.54 1.76 1.79 0.82 1.29

kinetic patterns.

Each reaction chamber is associated with

Semiconductor chip uses a filterless system that excludes excitation light based

Evanescent illumination at reaction chamber bottoms from nearby waveguide.

Functionalization of peptides at C-terminal lysines enables immobilization to the chip. Dve-labeled NAA recognizers and

Sequencing process proceeds in real time without fluidic exchange of reagents.

Recognition: 10s-100s of pulsing events

Fluorescence lifetime differentiates dye-

Aminopeptidases: Cleavage events stochastic at the single-trace level.

Signature summarizes the sequencing behavior of an Sensitive to AAs and variants such as

High-quality reads containing 4 or more RSs and 3 or more unique recognizers are eligible for alignment and aligned to a reference peptide sequence based on

Platinum reveals immunopeptides within a mixture of viral epitopes

MHC Class I Control Peptide Pool:

Mixture of 32 peptides from cytomegalovirus (CMV), Epstein-Barr virus (EBV), and influenza nucleoprotein (NP) virus. The pool consists of defined MHC class I-restricted T cell epitopes from these three viruses and can be used as a

- positive control to stimulate T cells.
- Bioinformatics indicates five sequenceable peptides on Platinum.

Viral peptide epitopes	So
RVLSFIKGTK	Influenza I
ILRGSVAHK	Influenza
RVRAYTYSK	EBVF
RLRAEAQVK	EBVF
IVTDFSVIK	EBV H

1. Elute, prepare, and sequence single peptide molecules on Platinum





CTELK.	5 T G G P I Y K*	R V L S F	1 K* +	о т к. 22 22	ILRGSVAHK*	
R V R A Y T Y S	SK ^a	A Q V K* *	. <mark>5 </mark>	P 5 G P L K*	50 A V F D R K ⁴	5 D A K* 60
I V T D F S V I	K ATIGT	A M Y K*	D Y C 1	NVLNK [®] 34	94 F L P F D K ^a 201	

We thank scientists Alexa Andrzejewski, Marla Charron, Cassie Simonides, Khanh D.Q. Nguyen, (Quantum-Si) and Professor Tanyeri Barak (Yale) for technical assistance.

References

Caron, Etienne, et al. "A case for a human immuno-peptidome project consortium." Immunity 47.2 (2017): 203-208

Kapoor, Saketh, et al. "Scaling up robust immunopeptidomics technologies for a global T cell surveillance digital network." Journal of Experimental Medicine 221.1 (2024). et al. "Real-time dynamic single-molecule protein sequencing on an integrated

emiconductor device." Science 378.6616 (2022): 186-192 kinner, Kenneth, David Kamber, and Brian Reed. "Detection of arginine posttranslational nodifications by single-molecule protein sequencing on the Quantum-Si platform." Cancer Research 83.7_Supplement (2023): 5301-5301.

Trademarks/Licensing

All trademarks are the property of Quantum-Si, Inc. or their respective owners. For specific trademark information, see: www.quantum-si.com/legal-disclaimer/

For research use only. Not for use in diagnostic procedures.